

1 mL Luer Lock Syringe

REF QIT011



3 mL Luer Lock Syringe

REF QIT012



20 mL Luer Lock Syringe

REF QIT013



INDICATIONS FOR USE

Qitexio Luer Lock Syringes are manual devices used to inject fluids into the body. They are also used for delivery of Lipiodol (Ethiodized Oil) Injection.





















PRECAUTIONS FOR USE

- This device must be used by qualified healthcare professionals only.
- This device is not made with natural rubber latex.
- Single-use device: this device is to be used once on one patient only.
- Do not re-sterilize. This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in bioincompatibility, infection or product failure risks to the patient.
- Before use, always check that the medical device and its packaging are undamaged.
- Before and during use, ensure all facility safety, hygiene, and asepsis protocol(s) are followed.
- For safe & proper use of the product, please carefully read & follow the full Prescribing Information (PI) for any drug product(s) used with the device. Ensure that the use of this device complies with the valid rules of good practice & current approved labeling of all drugs and devices used with this device.
- This device is intended for manual, hand-held use only. Do not connect to

active system(s).

- Connections are compatible with 6% (Luer) taper. The syringe barrel contains an ISO 80369-7 compliant fixed male connector, which is compatible with ISO 80369-7 compliant female luer hubs (ISO 80369-7: Small-bore connectors for liquids and gases in healthcare applications — Part 7: Connectors for intravascular or hypodermic applications).
- All Luer connections should be gently hand tightened without over tightening to ensure their security and to prevent damage to the device.
- Apply controlled axial pressure on the syringe's plunger to avoid uncontrolled injection and/or damage to the device.
- After use, the device must be handled and disposed of as hazardous biomedical waste in accordance with facility protocol(s) and any applicable federal, state, and local regulation(s).
- If the used medical device is to be returned to the manufacturer, ensure it is decontaminated first, according to the regulations in force for the device and the administered fluids.

DESCRIPTION OF SYMBOLS

	Consult instructions for use		Do not reuse		Do not re-sterilize
	Batch number		Sterilized using irradiation		Do not use if package open or damaged
	Manufacturer's product reference		Single sterile barrier system. Sterilized using irradiation		Keep dry
	Use-by date		Manufacturer		Fragile: handle with care
	Date of manufacture		Temperature limitation		Poses no known hazards in MR environments
	Verified for use with lipiodol® (Ethiodized Oil) Injection		Non-pyrogenic		Quantity (X) per package
	Prescription only		This way up		